


Data-Pak™ Battery




HEARTSINE™
TECHNOLOGIES

*A history of innovation...
A future of lifesaving technologies.*



Helping you... Deliver energy & store data

- 12 hours of monitoring
- 120 shocks at 200J
- 5 Year shelf life
- 2 Year stand-by
- Integrated flash memory for event & data recording



Long-lasting battery

The long-life Data-Pak™ battery provides over 12 hours of continuous monitoring or 120 full energy shocks, 5 year shelf-life or 2 year stand-by life. Each Data-Pak™ incorporates an integrated memory system that can record multiple patients, up to 12 hours of ECG including all event/incident recording and 60 minutes of audio.

Simple and Reliable Data Recovery

Each Samaritan® AED uses an Data-Pak™ Battery. This unique battery simultaneously powers the unit and records information (event resuscitation, patient ECG, logged events, timing and device status) which gives a complete patient record for physician oversight, reporting and Utstein data collection. Data transfer is as easy as removing the battery from the Samaritan and placing it into the Samaritan® Data® Recovery Unit (SDRU) which then downloads data to your PC. No more expensive data cards to lose, break or purchase.

Comprehensive Event Documentation

With the capacity to store up to 12 hours of ECG and 60 minutes of voice plus all associated text and logged events, the Saver® - Data-Pak™ - Samaritan® AED combination provides a complete documentary record of multiple resuscitations. Simple to install, easy to operate the Samaritan® data storage system provides one of the most comprehensive data review and storage system needed for first response and public access defibrillation programs.

Easily Upgrades one or more AED's

With Saver®, you can program an Data-Pak™ Battery with device configuration and upgrade data. These new configurations can then applied to the Samaritan® AED when the battery is installed into the unit. For example, you can use Saver® to set speaker volume, screen contrast, or time synchronization. The Data-Pak™'s reliable event documentation format prevents loss of critical information and helps you comply with AED reporting regulations.



1.866.HRT.SINE
(1.866.478.7463)

sales@heartsine.com
www.heartsine.com

Corporate Headquarters:
HeartSine Technologies, Inc.
940 Calle Amanecer - Suite E
San Clemente, CA 92673 - USA
Tel: 1.949.218.0092
Fax: 1.949.218.0093

European Operations:
HeartSine Technologies, Inc.
203 Airport Road West
Belfast, BT3 9ED - N. Ireland
Tel: +44 (0)28 9093 9400
Fax: +44 (0)28 9093 9401